



aycan myCD is a versatile digital media publisher solution for clinics, hospitals and imaging centres. The user can easily burn patient CDs just sending the DICOM images to **aycan myCD**. The software will automatically burn the single study CD, easy and fast. The authorized users can even burn multiple studies of a single patient or multiple studies of multiple patients in a CD/DVD with a simple and easy web based interface. The produced media is intended as a handout to patients only and not for diagnostic or therapeutic use.

Supporting the leading CD, DVD hardware (Epson, Rimage, Primera and Bravo), **aycan myCD** reduces the cost of film distribution, improves operational workflow and provides a highly professional product that improves patients services.

aycan myCD's workflow and burned content is fully compliant with IHE PDI requirements and produces a standard DICOM 3.0 part 10 formatted media. The studies can be sent from multiple DICOM sources including modalities, workstations or a PACS server.

The media contains DICOM compliant images in a DICOMDIR hierarchy. **aycan myCD** supports the printing of custom designed media labels, which can automatically include the patient and study information from the DICOM header of the images.



aycan myCD features:

- » Easy to install
- » Full automated process
- » Independent on media producer, tested with Epson, Rimage Primera, Bravo
- » Store and burn single patient study
- » Customizable media labels
- » Query and manage production process
- » Query and Retrieve for custom media
- » DICOMDIR directory hierarchy
- » Studies provided in HTML/JPEG format
- » Supports spanning of multiple media for large studies
- » Single patient, multiple study & multiple patient, multiple study media creation