

# C22WT

## 21.5" Multi-function Clinical Display



## 21.5" Multi-function Clinical Display

Fully compliant with the DICOM standard and with FHD widescreen design, C22WT can help doctors to access to patient information in a filmless environment while browsing different medical images. With excellent projective capacitive touch which supports finger, glove, stylus etc., C22WT can provide the perfect multi-touch function and handwritten signature system. That can be easily operated, and improve doctors' working efficiency and patient satisfaction. With the latest true flat design and high grade protective performance, C22WT can be easily cleaned and disinfected. Its durability is also increased.

**BEACON**

# C22WT

## 21.5" Multi-function Clinical Display

### Application Area

- Outpatient service
- Sickroom
- Intensive care
- Clinical department
- Community medical service



### Technical Specification

Model	C22WT
Display size	21.5"
Support color	16.7M
Resolution	1920 x 1080
Pixel Pitch	0.248 x 0.248 mm
Brightness	250 cd/m2 Typ
Contrast	1000:1 Typ.
Angle of view	178° Typ.
Response time	14 ms Typ.
Input signal interface	DisplayPort: 20pin DVI: DVI-D 24pin VGA: D-sub 15pin Audio: 3.5mm stereo jack Speaker: 5W*2
Size	523(W)x376(H)x179(D) mm
Weight	6.0kg
Standard mounting	VESA 100x100mm
Power supply	AC 100~240V/DC 12V
Power Consumption	40 max.
Power Saving	EPA, VESA
MTBF	> 50,000 Hours Support touch&electronic whiteboard



Professional protective glass



Built-in multiple display mode



Backlight stabilization system



Compliant with DICOM standard

**BEACON**

Shenzhen Beacon Display Technology Co., Ltd.

Room 201 Incubator Building, CASTD, High-tech

Tel: +86 755 2699 5355

South 1st Street, Nanshan District, Shenzhen 518057,

Fax: +86 755 2699 5755

China

www.beacon-display.cn

All specifications are subjected to change without prior notice.