


# C32SP+

3MP Color Medical Display



## 3MP Color Medical Display

C32SP+ is with 3MP (2048X1536) high resolution, high brightness, high contrast ratio and wide viewing angle, can be applied to various medical image equipment and PACS system.

Each monitor is strictly calibrated to compliant with DICOM standard, ensure the accurate display of medical images and details. The built-in backlight stabilization system can ensure the consistence performance in long-term use.

C32SP+ is with LED backlight, ultra-thin true flat design and AR protection glass, which makes it is anti-reflective, clean and disinfect easily.


**BEACON**

# C32SP+

## 3MP Color Medical Display

### Application Area

- CR/DR
- CT
- MRI
- PET/CT
- X-RAY
- DSA
- DSI
- 3D medical image
- Image fusion
- Nuclear medicine
- CAD



### Technical Specification

Model	C32SP+
Display size	21.3"
Support color	1.07 Billion
Resolution	2048x1536
Pixel Pitch	0.2115 x 0.2115mm
Brightness	800cd/m2 Typ.
Contrast	1400:1 Typ.
Angle of view	176° Typ.
Response time	20ms Typ.
Input signal interface	DVI-D: 24 pin DisplayPort: 20 pin
Size	369(W)x507~627(H)x235(D)mm
Weight	9.7 kg
Standard Mounting	VESA 100x100mm
Power supply	DC 12V
Power Consumption	60W max
Power Saving	EPA, VESA
MTBF	> 50,000 hours



Continuous quality assurance system(CQAS)



Backlight stabilization system



Ambient light adaptive system(ALAS)



QA software



Professional protective glass



Compliant with DICOM standard

**BEACON**

Shenzhen Beacon Display Technology Co., Ltd.

Room 201 Incubator Building, CASTD, High-tech  
South 1st Street, Nanshan District, Shenzhen 518057,  
China

Tel: +86 755 2699 5355

Fax: +86 755 2699 5755

[www.beacon-display.cn](http://www.beacon-display.cn)

All specifications are subjected to change without prior notice.