


VERT-X


High Performance Dedicated Chest X-ray Systems

LET'S START RIGHT NOW AND BE THE BEST IN YOUR BUSINESS

- Dedicated chest x-ray system with integrated moving grid bucky and DR detector designed for easy handling, control and high patient throughput.
- The professional DR acquisition software is user-friendly and can be quickly integrated into standard operating procedures, with integrated work-list function streamlines patient data entry.



Cost Effective DR Solution
System and function stability is extremely high since there are no mechanical parts like rollers, films etc. making this system virtually maintenance free.



Integrated MINI PACS
User friendliness to ensure quick and efficient workflow.



Improved Diagnosis
Due to optimal image quality and the option of computer assisted image manipulation.



INNOMED TOP-X 100LC HF X-RAY GENERATOR

High frequency X-ray generator that operates from conventional 220v 20 amp plug point.



Also compatible with any make of CR.