



TOP-X CR CLASSIC X-RAY SYSTEM

The TOP-X CR Classic floor mounted CR system with tube stand, floating table and wall-bucky stand with high frequency generator, normal speed 300k HU x-ray tube and collimator.



COST EFFECTIVE CR SOLUTION

System and function stability is extremely high since there are no mechanical parts like rollers, films etc. making this system virtually maintenance free.



INTEGRATED MINI PACS

Due to optimal image quality and the option of computer assisted image manipulation.



IMPROVED DIAGNOSIS

User friendliness to ensure quick and efficient workflow.



TOP-X 100LC HIGH FREQUENCY X-RAY GENERATOR



IMAGE X
Medical Imaging & Radiology Supplies

Get in touch:

- Tel: 011 869 6888
- www.imagex.co.za

S P E C I F I C A T I O N S

RADIOGRAPHIC TABLE	TUBE STAND
Size of Table top: 220cm by 81cm	Floor mounted tube stand
4-way floating table top with elevating movement (Vertical: 55-85cm)	Longitudinal travel: 195cm Vertical Travel: 162cm
X-way floating 55cm, Y-Way 12.5cm	Minimum focus to floor distance: 25cm & Maximum focus to floor distance: 187cm
Weight carrying capacity: 240kg	Electromagnetic brakes & Tube Horizontal rotation: 90 degrees
Built in bucky with 3-field AEC chamber and high-resolution grid	

WALL STAND	ROTATING ANODE X-RAY TUBE , MANUAL COLLIMATOR	
Floor mounted wall stand	Tube Voltage: 40-150kV	Anode Heat Capacity: 300kHU
Vertical Movement Range: 38-187cm	Focal Spots: 0.6/1.2mm (40/100 kW)	
Cassette Holder Bucky with 3-field AEC chamber and high-resolution grid	Speed of Rotation: 2800 Manual alignment to focus, on the collimator	
Cassette size: up to 17" x 17" (43x43cm)	Collimator Light source: LED	
	Extra flexible high voltage cables	

